EAI-EFAB43 – Rev.6 – 10/2020  
 ISO 5832-1; ISO 5832-2; ISO 5832-3;  
 ISO 5832-9; ISO 5834-2; ISO 13779-2

3S ORTHO  
 24, avenue Joannes Masset  
 FRANCE - 69009 Lyon  
 Tel : +33 4 37 24 07 45  
 Fax : +33 4 82 53 88 76  
 Email : [contact@3sortho.com](mailto:contact@3sortho.com)

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**English: ARAMIS/ INSTRUCTIONS FOR USE for the range ARAMIS total shoulder prosthesis. Important information for the use of 3S ORTHO Implants.**

**1-Generalities:**  
 Attention is drawn to the fact that the reliable function of implant is only guaranteed if the following basic rules are observed. It is important that the surgeon knows the surgical technique and that he takes account of these instructions for use. The patient selection is as important as the implants positioning.

**2-Description:**  
 The medical device ARAMIS is a total joint prosthesis designed to replace the glenohumeral joint of the shoulder in order to decrease pain and restore function and mobility in the shoulder articulation, comparatively to before the surgery.  
 Some osteosynthesis implants are made of many parts fitted together during the surgical operation. The implants are made of materials in accordance to European or International Standards in force. The composition of each implant together with the reference standards are indicated on the labels affixed on the individual packages. The surgeon should conduct tests if there is a risk that the patient might become allergic to the prosthetic materials. We strongly advise against using the 3S ORTHO implants with those of other companies.

**3-Materials:**  
 Humeral stem (cemented): titanium alloy (ISO 5832-3)  
 Cementless humeral stem: titanium alloy (ISO 5832-3) coated with microporous titanium (ISO 5832-2) and hydroxyapatite (ISO 13779-2)  
 Humeral head, humeral cup and glenoid head: stainless steel (ISO 5832-9)  
 Insert and glenoid component: polyethylene (ISO 5834-2)  
 Glenoid base: titanium alloy (ISO 5832-3) coated with microporous titanium (ISO 5832-2) and hydroxyapatite (ISO 13779-2)  
 Screw: titanium alloy (ISO 5832-3)

**4-Indications:**  
 The generally accepted arthroplasty indications which are common to the anatomical and reversed configuration are:

- Serious inflammatory arthropathy or osteoarthritis in case conservative therapy or alternative treatments failed or were considered unsuitable
- Arthropathy due to degenerative disease
- Acute trauma or traumaticum sequelae
- Previous failed joint replacement

The patient's joint must be anatomically and structurally suited to receive the selected implant(s). The choice of the prosthesis configuration (hemiarthroplasty, anatomical or reverse prosthesis) is the responsibility of the surgeon. However, priority must be given to the reversed configuration in the following cases:

- Eccentric omarthrosis
  - Rotator cuff tear
- 5-Contra-indications:**  
 • General or local infections or diseases • Major muscular, neuromuscular or vascular deficiencies in the affected limb • Serious deformity • Loss of bone or bone stock of poor quality which could create an unacceptable risk of implant instability • Any current illness that might affect the function of the implant • Risk of known or suspected allergies to the materials of implants mentioned at §3

**6-Factors which could lead to negative side effects, compromising the expected lifetime of the implant:**  
 • Serious osteoporosis • Considerable bone deformation • Systemic and metabolic disorders • Infectious diseases • Known drug addiction and/or drug medicine or alcohol dependency • High and/or repeated stress conditions due to excessive physical activity  
 The surgeon also must give the patient all information concerning the influence of these factors on the operation's success and give advice in order to minimize these risks. The patient must not only be mentally capable of realizing the significance of the precautions, but also of following them. The patient must consent to the operation and its associated risks.

**7-Possible negative side effects:**

- 7-1 Per operative  
 • Bone perforation or fracture
- 7-2 Postoperative  
 • Venous thrombosis and pulmonary embolism • Infections • Cardiovascular disorders • Haematomas and delayed wound healing • Allergic reaction • Loosening or migration of the implants can occur due to improper or incorrect implant fixation, trauma, malalignment, and/or excessive physical activity • Wear and/or deformation of articulating surfaces • Metallosis or osteolysis may be caused by wear debris
- 7-3 Specific  
 • Moving and/or loosening of one or all prosthetic components • Peri-articular adhesion and fibrosis • Implant breakage • Late deformation of one or all components due to defective fixation or improper positioning • Dislocations • Scapular notching  
 • Decreased range of motion • Unequal limb length due to improper positioning or selection of the implants • Secondary fractures • Capsulitis • scapular notching • Failure of consolidation • Incorrect positioning of the stem  
 The surgeon must inform the patient of possible interferences and/or complications associated with the implant during supplementary examinations (MRI)

**8-Precautions of use:**  
 Before implanting an implant, the surgeon must take into consideration the patient's overall clinical condition. The slightest efficiency or impossibility of any other treatments (non-operative and operative) have to be considered. He has also to take into account the previous surgeries performed on the patient. Protection against radiations. The surgeon will take all the necessary measures to ensure the protection towards the radiations caused by the per-operative scopic control of the bone fragments and the implants good positioning.

**9-Preoperative planning and preparation:**  
 Every implant is always subject to wear, which still remains a major problem awaiting for a solution. An initially stable implant may become loose in the course of time. Wear and loosening are two major causes that may drive to the re-operation of a failed implant. The surgeon must take into account the brochure and the operative manual which contain additional information about the implantation of the implants and the use of associated instruments. Their safety and functionality can only be warranted if the directions contained in the attached copy of the surgical procedure as well as the present warnings and instructions are duly considered and strictly observed. In case of any doubt, the surgeon must consult the supplier. Correct selection of the proper implant is of extreme importance and has to take into consideration the patient's weight, activity level and kind of bone default. It is recommended that surgeons use preoperative planning with the provided templates in order to determine the size of the prosthesis to implant.

**10-Implantation:**  
 The surgeon shall not use this prosthetic system with other components of other system and/or of other companies.  
 The surgeon has to respect the size compatibility between the different components; it is described in the following tables:

Glenoid	EAI G030	EAI G033	EAI G036	EAI G136
Humeral head				
EAI 4013	✓	✓	✓	X
EAI 4315	✓	✓	✓	X
EAI 4617	✓	✓	✓	X
EAI 4918	X	X	X	✓
EAI 4920	X	X	X	✓

Baseplate	EAI 0B28	EAI 0BHU	EAI 0BPS	EAI 0BRS	EAI 0BPL	EAI 0BRL
Glenosphère						
EAI SG38	✓	✓	✓	X	✓	X
EAI SG42	✓	✓	✓	X	✓	X
EAI SG35	X	X	X	✓	X	✓
EAI SC38	✓	✓	✓	✓	✓	X
EAI SC42	✓	✓	✓	X	✓	X
EAI SC35	X	X	X	✓	X	✓

It is recommended that surgeons use instruments which are specifically adapted to the 3S ORTHO components implantation.  
 The surgeon must ensure that:  
 • the implants are handled by people qualified for orthopaedic and traumatic operations  
 • the surgical implantation is realised under completely aseptic conditions.  
 The surgeon must verify the implants and ancillary's integrity before the operation. If the surgeon unintentionally damages the implant, packaging and/or instruments or notices any damage on them, they must not be used under no circumstances. Implants must not be reused: patient's contamination risks and mechanical properties non-guaranteed.

**11-Technical fixation:**  
 Proper fixation is paramount importance for the firm and long term seating of the implant. The following factors may give rise to complications:  
 • Excessive weakening of the cortical bone structure during the surgical preparation of the osseous bed • Selection of unsuitable implant size  
 • Excessive application of force during implant insertion or positioning, which may induce forced fractures or tears in the bone.  
 Humeral stems are indicated for either cemented or uncemented fixation applications. The product label specifies this aspect with two mentions: « cemented » or « cementless ». In case of cement use, it is essential to follow attentively the use conditions of the manufacturer. The cement manipulation can influence the implant fixation. The polyethylene glenoid components are indicated for cemented applications. Depending on the rasp stability, the surgeon must select the size and type appropriate for the patient.

**12-Post-operative monitoring:**  
 Clinical success depends on the validity of the indication, the selection of implants of the proper size, precise surgical procedure, the respect of the recommended conditions for use and good implant fixation. Post-operative monitoring periodic clinical and x-ray examinations in order to verify the artificial joint's performance is the surgeon's responsibility. If any complications and/or any dysfunction which might constitute a danger in regards to the patient's safety occur, another operation may become necessary. It is strongly recommended to the surgeon to inform the patient about advantages and disadvantages of the implant. The patient must be aware of the limitations of the implant and should also be warned that no

implant can be expected to withstand unsupported stress of full body weight. Non observance of these precautions may lead to serious consequences.

**13-Ancillary:**  
 ARAMIS shoulder prosthesis must be implanted with a dedicated instrument set. A surgical technique describing this specific instrument set as well as the "cleaning, sterilization" instruction (reference ISUP26) are enclosed at each delivery.

It is the user who verifies the status of the ancillary before surgery. If the user suspects or accidentally leads to a deterioration of the instrument set, these must be replaced. Non-sterile implants and instruments used for the implantation of this system must be sterilized before use. We recommended a steam sterilization method according to a protocol which is at least equivalent to what is required in the European Pharmacopoeia and is consistent with the national law in effect (for example: 134°C/138°C - 2 bars - 18 min). Instruments used for implantation must be carefully disinfected, cleaned and sterilized after each operation. Implants delivered sterile have been initially sterilized with a treatment validated conform to European or International Standards.

**14-Warnings:**  
 We recommend to check the state of the taper of the humeral stem in place, and to clean and dry it before the impaction of the component.  
 Store sterile implants in their original packages. The implants must be stored in a clean and dry area at room temperature. If either the implant or the package appear damaged, the implant should not be used.  
 The sterilization method is specified on the labels stuck on the individual packaging. Do not resterilize the implants.

**15-Information relating to magnetic resonance imaging (MRI):**  
 Generally, risks associated with a passive implant in a magnetic resonance (RM) environment include: heating, migration, and image artifacts at or near the implant site. Security, compatibility, warm-up and migration of ARAMIS system have not been evaluated in an RM environment. However, no complication has been reported to 3S ORTHO.  
 The surgeon must inform the patient of possible interferences and/or complications associated with the implant during supplementary examinations (MRI).  
 Other specific information relative to the different implants can be obtained at any time from 3S ORTHO, its subsidiaries or its distributors.

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