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English: ARAMIS-S INSTRUCTIONS FOR USE for the range ARAMIS-S total shoulder prosthesis. Important information for the use of 3S ORTHO Implants.

1-Generalités :

Attention is drawn to the fact that the reliable function of implant is only guaranteed if the following basic rules are observed. It is important that the surgeon knows the surgical technique and that he takes account of these instructions for use. The patient selection is as important as the implants positioning.

2-Description :

The medical device ARAMIS-S is a total joint prosthesis designed to replace the glenohumeral joint of the shoulder in order to decrease pain and restore function and mobility in the shoulder articulation, comparatively to before the surgery.

Some osteosynthesis implants are made of many parts fitted together during the surgical operation. The implants are made of materials in accordance to European or International Standards in force. The composition of each implant together with the reference standards are indicated on the labels affixed on the individual packages. The surgeon should conduct tests if there is a risk that the patient might become allergic to the prosthetic materials.

3-Matérials:

Humeral base (cementless) : titanium alloy (ISO 5832-3) coated with hydroxyapatite (ISO 13779-2)

Humeral centered head (cementless) : stainless steel (ISO 5832-9)

4-Indications:

Generally accepted indications for joint replacement include:

- Primary or secondary arthrosis of the glenohumeral joint where conservative therapy or alternative treatments have failed or are considered unsuitable
- Post-traumatic osteoarthritis
- Osteoarthritis due to instability
- Rheumatoid arthritis
- Humeral head avascular necrosis
- Previous surgery of glenohumeral joint that doesn't compromise the fixation of the new implant

The patient's joint must be anatomically and structurally suited to receive the selected implant(s). The choice of the prosthesis configuration (hemiarthroplasty or anatomical) is the responsibility of the surgeon.

5-Contra-indications:

- General or local deformities or diseases • Major muscular, neuromuscular or vascular deficiencies in the affected area • Serious deformity • Loss of bone or bone stock of poor quality which could create an unacceptable risk of implant instability • Any current illness that might affect the function of the implant • Risk of known or suspected allergies to the materials of implants mentioned at §3

6-Factors which could lead to negative side effects, compromising the expected lifetime of the implant:

• Serious osteoporosis • Considerable bone deformation • Systemic and metabolic disorders • Infectious diseases

• Known drug addiction and/or drug medicine or alcohol dependency

• High and/or repeated stress conditions due to excessive physical activity

The surgeon also must give the patient all information concerning the influence of these factors on the operation's success and give advice in order to minimize these risks. The patient must not only be mentally capable of realizing the significance of the precautions, but also of following them. The patient must consent to the operation

and its associated risks.

7-Possible negative side effects:

7-1 Per operative

• Bone perforation or fracture

7-2 Postoperative

• Venous thrombosis and pulmonary embolism • Infections • Cardiovascular disorders • Haematomas and delayed wound healing • Allergic reaction • Loosening or migration of the implants can occur due to improper or incorrect implant fixation, trauma, malalignment, and/or excessive physical activity • Wear and/or deformation of articulating surfaces • Metallosis or osteolysis may be caused by wear debris

7-3 Specific

• Moving and/or loosening of one or all prosthetic components • Peri-articular adhesion and fibrosis • Implant breakage • Late deformation of one or all components due to defective fixation or improper positioning • Dislocations • Scapular notching

• Decreased range of motion • Unequal limb length due to improper positioning or selection of the implants • Secondary fractures • Capsulitis • Scapular notching • Failure of consolidation • Incorrect positioning of the stem

The surgeon must inform the patient of possible interferences