


ARAMIS-S






TOTAL SHOULDER PROSTHESIS

Performance

Security



- Less constraints on the diaphyseal zone of the humerus
- Maximal bone conservation
- Anatomic adaptability
- Simple procedure
- Stemless and Cementless (HA coating for Bone Integration)


-  Primary or secondary arthritis of the glenohumeral joint where conservative therapy or alternative treatments have failed or are considered unsuitable
-  Post-traumatic with or without vicious cal
-  Rheumatoid arthritis
-  Humeral head avascular necrosis
-  Previous surgery of glenohumeral joint that does not compromise the fixation of the new implant

3S

O R T H 

ARAMIS-S


TOTAL SHOULDER PROSTHESIS



Humeral head resection




Centering pin



Humeral preparation



Thread with a tap



Placement of the implant



**Simple et Precise
Instrumentation**



Reference

ARAMIS-S system

- ESI2122** Humeral base Ø21/L22mm
Cementless
- ESI2124** Humeral base Ø21/L24mm
Cementless
- ESI4014** Humeral centered head Ø40
H14mm Cementless
- ESI4315** Humeral centered head Ø43
H15mm Cementless
- ESI4617** Humeral centered head Ø46
H17mm Cementless
- ESI4918** Humeral centered head Ø49
H18mm Cementless
- ESI4920** Humeral centered head Ø49
H20mm Cementless



3S
ORTHO

Class III medical device // Indication: joint arthroplasty of the glenohumeral joint



Consult the surgical technique before use



Read the instructions for use carefully



www.3sortho.com



Les passerelles

24 avenue Joannès Masset

69009 LYON, FRANCE

Phone : +33 4.37.24.07.45

Fax : +33 4.72.74.90.41

contact@3sortho.com